



スキャンラン

ティップ・ガード (薬事非該当品)

手術用器具及び内視鏡先端保護キャップ

- カラーコート化された特殊ポリマー製キャップです。
- キャップしたまま蒸気やガス滅菌可能です。
- X線不透過です。
- カッチューブより有効・能率的です。
- 器具用 1袋100ヶ入
- 内視鏡用 1袋30ヶ入

CATALOG NUMBERS COLORED TIP-GUARD™ INSTRUMENT PROTECTORS



1001-89

WIDTH 1/8" 3.2mm
LENGTH 1" 25.4mm
HEIGHT 1" 25.4mm



1001-88

1/16" 1.6mm
1" 25.4mm
5/8" 15.9mm



1001-87

1/16" 1.6mm
1" 25.4mm
3/8" 9.5mm



1001-91

13/32" 10.3mm
3/4" 19mm



1001-90

3/16" 4.8mm
1" 25.4mm



1001-85

1/8" 3.2mm
1" 25.4mm




1001-86

7/64" 2.8mm
3/4" 19mm



1001-93

5/64" 2mm
3/4" 19mm



1001-92

1/16" 1.6mm
3/4" 19mm

TINTED TIP-GUARD™ INSTRUMENT PROTECTORS



1001-305

WIDTH 1/8" 3.2mm
LENGTH 1" 25.4mm
HEIGHT 1" 25.4mm



1001-304

1/16" 1.6mm
1" 25.4mm
5/8" 15.9mm



1001-303

1/16" 1.6mm
1" 25.4mm
3/8" 9.5mm



1001-307

13/32" 10.3mm
3/4" 19mm



1001-306

3/16" 4.8mm
1" 25.4mm



1001-301

1/8" 3.2mm
1" 25.4mm




1001-302

7/64" 2.8mm
3/4" 19mm



1001-309

5/64" 2mm
3/4" 19mm



1001-308

1/16" 1.6mm
3/4" 19mm

CLEAR TIP-GUARD™ INSTRUMENT PROTECTORS



1001-205

WIDTH 1/8" 3.2mm
LENGTH 1" 25.4mm
HEIGHT 1" 25.4mm



1001-204

1/16" 1.6mm
1" 25.4mm
5/8" 15.9mm



1001-203

1/16" 1.6mm
1" 25.4mm
3/8" 9.5mm



1001-207

13/32" 10.3mm
3/4" 19mm



1001-206

3/16" 4.8mm
1" 25.4mm



1001-201

1/8" 3.2mm
1" 25.4mm



1001-202

7/64" 2.8mm
3/4" 19mm



1001-209

5/64" 2mm
3/4" 19mm



1001-208

1/16" 1.6mm
3/4" 19mm

製造販売元



日本エー・シー・ピー株式会社

〒113-0033 東京都文京区本郷2-27-3 菊屋ビル2階
電話 (03) 3814-2345 (代) ファックス (03) 3814-7484